



Worth a New Look

RELIANCE LOW VOLTAGE ENT TREATMENT CABINETS

ENT TREATMENT CABINETS

Model 310



Model 310 Narrow
in Shadow Laminate

Overall Height: 48"
Height to Work Surface: 35"
Width: 25"
Depth: 20 7/16"
Shipping Weight: 220 lbs.

Top Drawer:
8 13/16" W x 10" L x 1 3/4" D
Middle Drawers (2):
20 1/2" W x 12 1/2" L x 3" D
Bottom Drawers:
20 1/2" W x 12 1/2" L x 4" D

- Economical, narrow footprint, laminated treatment cabinet
- Spacious work surface
- 1200cc disposable suction container
- Two Welch Allyn instruments
- Three pieces of glassware set
- Vacuum and pressure pumps and hoses
- Black laminate top standard

Model 410



Model 410
in Natural Pear Laminate

Overall Height: 48"
Height to Work Surface: 35"
Width: 33"
Depth: 20 7/16"
Shipping Weight: 245 lbs.

Top Drawer:
12 13/16" W x 10" L x 1 3/4" D
Middle Drawers (2):
28 1/2" W x 12 1/2" L x 3" D
Bottom Drawers:
28 1/2" W x 12 1/2" L x 4" D

- Traditional laminated treatment cabinet – choose your own laminate color
- Easy accessibility to instruments
- 1200cc disposable suction container
- Generous drawer space
- Two Welch Allyn instruments
- Vacuum and pressure pumps and hoses
- Your choice of three pieces of glassware
- Black laminate top standard

Model 510



Model 510
in Stainless Steel

Overall Height: 48"
Height to Work Surface: 35"
Width: 33"
Depth: 20 7/16"
Shipping Weight: 275 lbs.

Top Drawer:
12 13/16" W x 10" L x 1 3/4" D
Middle Drawers (2):
28 1/2" W x 12 1/2" L x 3" D
Bottom Drawers:
28 1/2" W x 12 1/2" L x 4" D

- Stainless steel treatment cabinet
- Spacious work surface
- 1200cc disposable suction container
- Two Welch Allyn instruments
- Vacuum and pressure pumps and hoses
- Your choice of three pieces of glassware
- Black laminate top standard

Model 610



Model 610 Tall
in North Sea Laminate

Overall Height: 58"
Height to Work Surface: 45"
Width: 33"
Depth: 20 7/16"
Shipping Weight: 300 lbs.

Top Drawer:
12 13/16" W x 10" L x 1 3/4" D
Middle Drawers (2):
28 1/2" W x 12 1/2" L x 3" D
Bottom Drawers:
28 1/2" W x 12 1/2" L x 4" D
Extra Drawer:
28 1/2" W x 12 1/2" L x 4" D

- Tall traditional treatment cabinet with extensive storage
- Easy accessibility to instruments
- Generous drawer space
- 1200cc disposable suction container
- Two Welch Allyn instruments
- Vacuum and pressure pumps and hoses
- Your choice of three pieces of glassware
- Black laminate top standard

STANDARD AND OPTIONAL FEATURES

Standard components of all basic models of Reliance ENT treatment cabinets include:



Heavy duty Gas pressure & suction pumps

- Feature a 5 micron air filter and air pressure regulator that eliminates particulates



Elevated work tower

- Includes 3 pieces of glassware and 2 stainless steel canisters with lids
- Allows easy access to instruments and storage for dispensing gloves and tissues



Low voltage, backlit membrane control switches

- Meets IEC 60601-1 safety requirements
- Available in right hand or left hand configuration



Detachable power and footswitch cords

- Conveniently located in the back of the cabinet
- Recessed panel and cords save valuable floor space



Welch Allyn instruments

- Your choice of two rechargeable or corded instruments
- Can be mounted on either the left or right side of the cabinet



Full-extension drawers

- Removable inserts for easy cleaning
- Allows for versatile arrangement of instruments



Conveniently located hose clips

- Mount to either side of the cabinet providing a custom configuration for any practice



Sturdy, lockable 3-inch casters

- Allow for greater mobility
- Provide improved cabinet stability



Antimicrobial, curved drawer and door handles

- Standard for improved safety and style

Optional components of all basic models of Reliance ENT treatment cabinets include:



Three-port, LED light source

- Built-in option frees up valuable counter space
- Accepts five popular standard fittings
- Side mount only on the model 310



Premier bead warmer

- Can be recessed into the cabinet top
- Provides added safety and convenience



Corded, detachable footswitch

- Allow hands-free operation of suction and pressure pumps
- Can be added to any model at any time



Locking drawers

- Additional measure of safety and compliance
- Keyed and keyless entry are both available



Stainless steel debris tray and lid

- Can be mounted to the side of any cabinet model
- Tray and lid are easily removed and are autoclave safe

Add a concealed drawer warmer to keep mirrors fog-free. (Not pictured.)

LAMINATE AND SURFACE OPTIONS

Laminate Surface Options

Cabinet sides, backsplash, doors and drawers.



Desert Zephyr



Dove Grey



Natural Pear



North Sea



Satin Stainless



Shadow




Stainless Steel



White Nebula


Counter Surface Options

Standard Surface



Black Laminate

Solid Surface Options (Call for pricing)



Cloudy



Hazelnut



Cocoa Vine



Eclipse

Electrical Requirements: operates at 120V/60Hz and 220-240V 50-60Hz

The cabinets conform to the following standards:

ANSI/AAMI IEC60601-1: 2005

CAN/CSA CSA-C22.2 No. 60601-1-08