



# Worth a New Look

## RELIANCE ENT PRODUCT BROCHURE

## Treatment / Procedure Chairs



### 980 Full Power Procedure Chair

- Round base increases accessibility to patient
- Stylish removable cushions
- Detached footswitch controls allow for utilizing as both examination and procedure chair
- Programmable side switches
- Choice of headrests
- Low model assists ADA compliance



### 7000 Procedures Chair

- Detached full footswitch controls allows for utilizing as both examination and procedure chair
- Fully powered base and chair top recline
- Programmable side switches
- Choice of headrests
- Low model assists ADA compliance

### FX 920 Fully Powered Tilt Chair

- Patented motorized tilt function via membrane switches located on both sides of chair back or through the corded footswitch
- Efficient memory recline positioning
- Tilt design ideal for tight spaces
- Also available in manual recline
- Low model assists ADA compliance



### 710 Minor Procedures Chair/ Table with Pillow Headrest

- Reclines 180 degrees for minor procedures
- Optional foldaway headrest
- Memory positioning and safety switch



### 6200 Examination Chair

- Power hydraulic base and an assisted manual recline
- Reclines to full flat position
- Corded footswitch controls raising and lowering of base
- Membrane switches on both sides of chair control raise/lower/auto return function
- Low model assists ADA compliance

## Low Voltage ENT Treatment Cabinets

### Model 310



Model 310 Narrow  
in Shadow Laminate

Overall Height: 48"  
Height to Work Surface: 35"  
Width: 25"  
Depth: 20 7/16"  
Shipping Weight: 220 lbs.

Top Drawer:  
8 13/16" W x 10" L x 1 3/4" D  
Middle Drawers (2):  
20 1/2" W x 12 1/2" L x 3" D  
Bottom Drawers:  
20 1/2" W x 12 1/2" L x 4" D

- Economical, narrow footprint, laminated treatment cabinet
- Spacious work surface
- 1200cc disposable suction container
- Two Welch Allyn instruments
- Three pieces of glassware set
- Vacuum and pressure pumps and hoses
- Black laminate top standard

### Model 410



Model 410  
in Natural Pear Laminate

Overall Height: 48"  
Height to Work Surface: 35"  
Width: 33"  
Depth: 20 7/16"  
Shipping Weight: 245 lbs.

Top Drawer:  
12 13/16" W x 10" L x 1 3/4" D  
Middle Drawers (2):  
28 1/2" W x 12 1/2" L x 3" D  
Bottom Drawers:  
28 1/2" W x 12 1/2" L x 4" D

- Traditional laminated treatment cabinet – choose your own laminate color
- Easy accessibility to instruments
- 1200cc disposable suction container
- Generous drawer space
- Two Welch Allyn instruments
- Vacuum and pressure pumps and hoses
- Your choice of three pieces of glassware
- Black laminate top standard

### Model 510



Model 510  
in Stainless Steel

Overall Height: 48"  
Height to Work Surface: 35"  
Width: 33 "  
Depth: 20 7/16"  
Shipping Weight: 275 lbs.

Top Drawer:  
12 13/16" W x 10" L x 1 3/4" D  
Middle Drawers (2):  
28 1/2" W x 12 1/2" L x 3" D  
Bottom Drawers:  
28 1/2" W x 12 1/2" L x 4" D

- Stainless steel treatment cabinet
- Spacious work surface
- 1200cc disposable suction container
- Two Welch Allyn instruments
- Vacuum and pressure pumps and hoses
- Your choice of three pieces of glassware
- Black laminate top standard

### Model 610



Model 610 Tall  
in North Sea Laminate

Overall Height: 58"  
Height to Work Surface: 45"  
Width: 33"  
Depth: 20 7/16"  
Shipping Weight: 300 lbs.

Top Drawer:  
12 13/16" W x 10" L x 1 3/4" D  
Middle Drawers (2):  
28 1/2" W x 12 1/2" L x 3" D  
Bottom Drawers:  
28 1/2" W x 12 1/2" L x 4" D  
Extra Drawer:  
28 1/2" W x 12 1/2" L x 4" D

- Tall traditional treatment cabinet with extensive storage
- Easy accessibility to instruments
- Generous drawer space
- 1200cc disposable suction container
- Two Welch Allyn instruments
- Vacuum and pressure pumps and hoses
- Your choice of three pieces of glassware
- Black laminate top standard

## Examination & Surgical Stools

### 4200 Series Stools

Reliance's Economy Stools feature:

- Five legs to distribute weight evenly and ensure stability
- Height adjusts by pneumatic raising and lowering with lever
- 23" diameter base



4240  
basic

4246  
with backrest

### 5300 Pneumatic Stools

Reliance's top-of-the-line Pneumatic Examination Stools feature:

- Deep polyfoam seat offers support and comfort during lengthy procedures
- Five legs to distribute weight evenly and ensure stability
- 5300 Series height adjusts by using the control ring located under the seat
- Upholstered three-way, wrap-around arm for abdominal, back or arm support (5356 only)



5340  
basic

5346  
with backrest

5348  
waterfall cushion

5356  
wrap-around arm

### 500 Hydraulic Stools

Reliance's Hydraulic Surgical Stools Feature:

- Hydraulic raising and lowering foot pedal
- Deep polyfoam seat offers support and comfort during lengthy procedures
- Five legs to distribute weight evenly and ensure stability
- Floor lock standard on model 558, optional on all other hydraulic stools



540  
basic

546  
round seat

548  
waterfall cushion



556  
wrap-around arm

558  
individual arm rests

## ENT Light

### Halogen Examination Light

- Chair mounted examination light attaches to Reliance models 980, 710, 7000 and 6200
- Features a 30" swivel head gooseneck
- Light on/off switch is located on head of the examination light
- 20 watt, high intensity lamp

