



LASER RESISTANT ENDOTRACHEAL TUBE



Tenax™





# Tenax™

## LASER RESISTANT ENDOTRACHEAL TUBE

Tenax™ has an aluminum wrapped, reinforced spiral shaft that runs the length of the tube. Latex-free, Tenax is encased in a smooth sheath of flexible silicone, minimizing damage to vocal cords and related anatomy.

The Tenax™ DualCuff™ offers an extra layer of protection and is preloaded with non-staining blue dye.

### SPECIFICATIONS

PRODUCT NUMBER	INNER DIAMETER	OUTER DIAMETER	LENGTH
TG0050-S	5.0mm	7.9mm	34cm
TG0055-S	5.5mm	8.6mm	34cm
TG0060-S	6.0mm	9.3mm	34cm
TG0065-S	6.5mm	10.0mm	34cm
TG0070-S	7.0mm	10.6mm	34cm
TG0075-S	7.5mm	11.0mm	34cm

Each Tenax™ comes with stylet

