




THE FUTURE OF INTUBATION

**EMS, EMERGENCY MEDICINE, ANESTHESIA,
COMBAT MEDICS, ENT**

“Improve your view and your access to the trachea.”

www.adroitsurgical.com

About The Company



ADROIT SURGICAL

Adroit Surgical, LLC, is a privately funded medical device company. Adroit has developed a patented method for securing the difficult airway in emergency situations that is an alternative to a video laryngoscopy and in some cases cricothyrotomy. The first two introductory products are disrupting the current standard of care, making intubation faster, safer and less expensive. The Vie Scope® offers Paramedics and Doctors the unique ability to handle all adult intubations using only one device. The direct line of sight view for the user will lower patient complications as well as saving time, money and most importantly lives.

Why Choose Us?

At Adroit Surgical, we understand the critical importance of securing the airway swiftly and efficiently. As a leading medical device company, we specialize in providing innovative intubation devices that prioritize patient safety and improve procedural efficiency.



ENHANCED LINE OF SIGHT VIEW



EASY TO USE



PATENTED LIFESAVING TECHNOLOGY



Our Devices



The Vie Scope is a self-contained, battery-powered, disposable scope that takes advantage of a closed circular tube with a beveled end to visualize the vocal cords better than all other methods using the naked eye. The ring of LED lights is transmitted through the side wall of the tube down the entire length of the tube to give the user maximal illumination of the target tissue with minimal chance of light obstruction by secretions or blood.



Bébé Vie Scope™, developed by Adroit Surgical, is an innovative medical device designed specifically for intubating toddlers up to 10 years of age. This self-contained and battery-powered laryngoscope offers a convenient and efficient solution for healthcare professionals working in pediatric care.



The technique for intubation can be made safer, faster, and more efficient with the use of our Voir Bougie® (patent pending). The Voir Bougie® is unique in that clearly marked bands permit the user to note the correct depth upon insertion to avoid lung or tracheal injury. This allows the process of intubation to occur safely when passing the endotracheal tube over the Voir Bougie® during intubation.


TRUSTED AND RELIABLE TRIO TO
SAVE LIVES



Vie Scope®




bébé Vie Scope™



Voir Bougie®

Contact Us

Now!



@adroitsurgical
www.adroitsurgical.com

1-877-7ADROIT

Adroit Surgical LLC- 1000I Broadway
Extension Suite A, Oklahoma City, OK
73114, USA