



Preserving What Matters

SUBDERMAL NEEDLE ELECTRODES

Reliable Response Acquisition
For Dynamic Electromyography
Recording During Surgery

- ✓ Same color coding from patient to display
- ✓ Packaged in sets to improve efficiency
- ✓ Suitable for use with multi-channel nerve monitors
- ✓ 14 mm subdermal needle
- ✓ 1.5 m twisted pair cables, which connect to the pre-amplifier pod for ease of use
- ✓ Single patient use to reduce cross-contamination risk

Neurosign® disposable subdermal needle electrodes are designed to acquire EMC signals from muscles. When used with the Neurosign® V4 Intraoperative Nerve Monitor, the colors carry right through from connection to the on-screen channels.



Manufacturer:

Technomed
Amerikalaan 71
6199 AE Maastricht-Airport
The Netherlands
Tel: +31 43 408 68 68

Magstim Inc., dba Neurosign

17200 Medina Road
Ste. 800
Plymouth, MN 55447
844.MAGSTIM (844.624.7846)
info@neurosign.com



SUBDERMAL NEEDLE ELECTRODES

<p>2-Channel Electrode Set PN 4469-00</p> <p>The 2-channel electrode set contains all required needle electrodes for 2-channel EMG monitoring, including a reference needle electrode (green) and a return needle electrode (white) for use with stimulating probes.</p>	<p>This electrode set includes:</p> <ul style="list-style-type: none"> • Blue electrodes, twisted pair (Ch1) • Red electrodes, twisted pair (Ch2) • Green electrode, single (Reference) • White electrode, single (Return)
<p>4-Channel Electrode Set PN 4470-00</p> <p>The 4-channel electrode set contains all required needle electrodes for 4-channel EMG monitoring, including a reference needle electrode (green) and a return needle electrode (white) for use with stimulating probes.</p>	<p>This electrode set includes:</p> <ul style="list-style-type: none"> • Blue electrodes, twisted pair (Ch1) • Red electrodes, twisted pair (Ch2) • Violet electrodes, twisted pair (Ch3) • Orange electrodes, twisted pair (Ch4) • Green electrode, single (Reference) • White electrode, single (Return)
<p>8-Channel Electrode Set PN 4472-00</p> <p>The 8-channel electrode set contains all required needle electrodes for 8-channel EMG monitoring, including a reference needle electrode (green) and a return needle electrode (white) for use with stimulating probes.</p>	<p>This electrode set includes:</p> <ul style="list-style-type: none"> • Blue electrodes, twisted pair (Ch1) • Red electrodes, twisted pair (Ch2) • Violet electrodes, twisted pair (Ch3) • Orange electrodes, twisted pair (Ch4) • Grey electrodes, twisted pair (Ch5) • Yellow electrodes, twisted pair (Ch6) • Brown electrodes, twisted pair (Ch7) • Beige electrodes, twisted pair (Ch8) • Green electrode, single (Reference) • White electrode, single (Return)

When used with the Neurosign® V4 Intraoperative Nerve Monitor, the color coded needle electrodes correspond to both the pre-amplifier pod and the EMG channels displayed on-screen.