

STIMULATING INSTRUMENTS



STIMULATE WHILE YOU OPERATE TO IMPROVE PATIENT OUTCOMES

- ✓ Ready to use
- ✓ 100% electrical and insulation reliability
- ✓ No insulation degradation when compared to reusable instruments
- ✓ Connection cable included
- ✓ Single patient use, supporting patient safety by reducing cross-contamination risk
- ✓ Cost-effective intraoperative nerve monitoring solutions

Neurosign® disposable stimulating instruments are designed to manipulate, lift, and dissect tissue while simultaneously stimulating during procedures such as spinal nerve roots or neuro-otology surgery. The built-in stimulation capabilities provides ongoing feedback of nerve proximity, improving efficiency and aiding in patient safety by mitigating the need to switch between tools.



REDUCING SURGERY TIME AND SUPPORTING PATIENT SAFETY AND QUALITY OF CARE

- **COMPOSITION** — machined from quality stainless steel
- **DESIGN** — sized for precision needs
- **SINGLE PATIENT USE INSTRUMENT** — individually packaged and sterile
- **ECONOMICAL** — defined cost per patient
- **REDUCES RISK** — Single patient use
- **QUALITY** — Smooth working surfaces when starting procedures.

Compatible with Neurosign® V4, N100, N400 and N800.

3MM RING DISSECTOR

The Ring Dissector is often used to remove a tumor inside tissue.

PN 5877-00

individually packed, sterile, single use, including a connection cable



CURVED NEEDLE

The Curved Needle is used to reach around tissue or separating nerves from a tumor.

PN 5876-00

individually packed, sterile, single use, including a connection cable



ELEVATOR

The Elevator is designed to elevate nerves from a tumor.

PN 5875-00

individually packed, sterile, single use, including a connection cable



Neurosign® is a registered trademark of The Magstim Company Limited. All specifications are subject to change. The information contained within this document is accurate at the time of publication.

© 2022 Technomed
N-2022-005.01_US

Legal Manufacturer:
Technomed Europe
Amerikalaan 71
6199 AE Maastricht-Airport
The Netherlands

Tel: +31 43 408 68 68
info@neurosign.com

neurosign.com

U.S. Headquarters:
Magstim Inc., dba Neurosign
17200 Medina Road
Ste. 800
Plymouth, MN 55447

Tel: 844.624.7846

Your NEUROSIGN® Business Partner: