

HGN KIT 1

PART NUMBER: NM-KIT 1



TWISTED PAIR ELECTRODES

For use with the Neurosign 400 and 800. These electrodes ensure excellent recordings whilst the twisted cable minimizes noise.

SINGLE NEEDLE ELECTRODE

For use as the ground/reference electrode with the Neurosign 100, 400 and 800 and as a return electrode for monopolar probes

BIPOLAR PROBE

The Bipolar Probe has tips separated by 1 mm and is ideal for use when stimulating through tissue and bone; for example during thyroid surgery.

PN 1699-00-TE | 20 mm x 0.3 mm
needle electrode.
1 m wire. Box of 20



PN 1705-00-TE | 20 mm x 0.3 mm
needle electrode
1 m wire. box of 20



PN 3601-00-TE | 100 mm shaft
length, 0.75 mm diameter,
2 mm exposed tip, 1mm tip
separation, box of 10

